

Paraffin oven UNpa with TwinDISPLAY AtmoCONTROL software


Model sizes: 30 / 55 / 75 / 110 / 160 +20 °C to +80 °C

PARAFFIN OVEN UNpa Five model sizes, five times high-precision temperature control of the embedding medium paraffin in science and research. The range of functions and thermal safety of paraffin ovens UNpa are designed specifically for absolutely reliable sample preparation in the laboratory. The benefits for the user: an optimal cost/benefit ratio for an appliance that guarantees, for many years, precise and even temperature control for embedding media without any loss in quality whatsoever.



Safe warming of paraffin

Thanks to its high capillarity, liquid paraffin is an ideal embedding medium. This property, however, may lead to oily residue in tiny cavities. For this reason, the interior chamber of paraffin ovens UNpa is designed almost gas tight. There is definitely no danger of ignition of residue or damage to mechanical and electronic components.



Absolutely uniform temperature distribution

Due to the almost gas tight chamber, no outside air is exchanged. Therefore, the advantages of the uniform temperature distribution by the large surface all-round heating system applied in Memmert heating ovens come fully into play. Also without forced convection, the perfect interaction of the control system and heating unit ensures unparalleled temperature homogeneity and stability.



Air flow with natural convection



PARAFFIN OVENS UNpa

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

Standard units are safety-approved and bear the test marks: C € 🚉 ເພື່ອເຂື້ອ







Interior:

Stainless steel, material 1.4301 (ASTM 304), with all-round deep-drawn ribs to integrate the largearea heating with ceramic-metal sheath, nearly gastight

Housing:

Textured stainless steel, rear zinc-plated steel, intuitively operated TwinDISPLAY (TFT colour display) with touchscreen, fully insulated stainless steel door


Mains cable with plug (German type) Connection:


Installation: 4 feet

Interfaces:









Model sizes/Descrip	otion		30	55	75	110	160	
Stainless steel	Volume	approx. l	32	53	74	108	161	
interior	Width	(A) mm	400			560		
	Height	(B) mm	320	400	560	480	720	
	Depth	(C) mm	250	33	30	4	00	
	Max. number of grids/shelves	number	3	4	6	5	8	
	Max. loading per grid/shelf	kg	20					
	Max. loading of chamber	kg	60	80	120	175	210	
	Max. loading per slide-in drip tray	kg	1.5				3	
	Max. loading per bottom drip tray	kg		1.5		3		
Textured stainless	Width	(D) mm		585		7	45	
steel exterior	Height	(E) mm	704	784	944	864	1104	
	Depth (without door handle, door handle +56 mm)	(F) mm	434	51	14	5	84	
Standard equipment	Stainless steel grids, electropolished	number	1 2					
	Standard works calibration certificate (measuring point chamber center)	°C	+80					
Temperature	Working temperature range	°C	at least 5 above ambient temperature to +80					
·	Setting temperature range	°C	+20 to +80					
	Setting accuracy	°C	0.1					
Further data	Electrical load at 230 V, 50/60 Hz	approx. W	1600	2000	2500	2800	3200	
	Electrical load at 115 V, 50/60 Hz	approx. W	1600					
Packing data	Net weight	approx. kg	45	55	66	75	96	
·	Gross weight (packed in carton)	approx. kg	61	74	85	100	122	
	Width	approx. mm	660 730			830		
	Height	approx. mm	890	950	1130	1050	1300	
	Depth	approx. mm	650 670		800			
Order No. Paraffin Ovens			UN30pa	UN55pa	UN75pa	UN110pa	UN160pa	

Options		30	55	75	110		160	
/oltage 115 V, 50/60 Hz				X2				
Full-sight glass door (4-layer insulating glass)				BO				
Entry port, 23 mm clear diameter, for	left centre/centre			F0				
ntroducing connections at the side, gastight, can be closed by flap and silicone stopper,	left centre/top			F1				
an be closed by hap and sincone stopper, tandard positions - technical clarification	right centre/centre			F2				
equired	right centre/top			F3				
Entry port, 23 mm clear diameter, gas tight,	left							
an be closed by flap and silicone stopper, in	right			F5				
pecial positions (please state location) - echnical clarification required	rear			F6				
<u>'</u>								
Entry port (silicone), 40 mm clear diameter, gas tight, can be closed bositions at the back (please state location) - technical clarification re	quired			F7				
. 20 4)	erature controller, actual value			V3				
flexibly in chan	of a Pt100 sensor positioned ober for external temperature ing (max. 3) - price per sensor	V6						
Gas inlet/outlet through 2 tubes with ball valves - technical clarificati vith other options	on required for combination	K3						
Norks calibration certificate for 3 temperatures: +37, +52, +70 °C				D00126				
Norks calibration certificate for one (freely selectable) temperature va specification	alue according to customer			D00109				
Door with lock and key (safety lock)				В6				
Door hinged on the left				В8				
otential-free contact for combination error message (e.g. supply fail	ure, sensor fault, fuse)			H6				
Potential-free contact (24 V/2 A) with socket, or signal generation, controlled by programme segment, for free-selectable unctions to be activated (e.g. activation of audible and visual signals, exhaust motors, ans, stirrers, etc.)	2 contacts			H72				
Process-dependent programmable door lock				D4				
Door-open-recognition		V5						
Flexible Pt100 temperature sensor, positioned flexibly in chamber or measurement (up to 1 additional sensor is possible). The measured to ndicated on the display, recorded in the integral data store, and can AtmoCONTROL software	emperature will be			Н8				
MobileALERT, notification by SMS in case of any error or alarm of the	device (requires option H6)			C3				
Castor frame (2-part), height 140 mm				R9				
Assessarios			າ			110		

Castor name (2 part), neight 140 min		II.J				
Accessories	30	55	75	110	160	
Stainless steel grid, electropolished	E28884	E28884 E20164 E2		E20	20165	
Perforated stainless steel shelf	B29727 B03916		B00	B00325		
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution)	E02070 E02072			E02	E02073	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution)	B04356 B04358 B0			B04	359	
Wall bracket for wall mounting	B29755	B29756	B29757	B29758	B29759	
Guarantee extension by 1 year		GA1Q5				
USB-Ethernet adapter			E06192			
Ethernet connection cable 5 m for computer interface			E06189			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number	B33170					
Set of height adjustable feet (4 pcs)	B29768					
Stacking set (4 pcs) for stacking of appliances of same size	B29744				-	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29732	B29734	B29736	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729	B29731	B29733	B29735	B29737	
Subframe, adjustable in height (size 30 to 75: height 600 mm, size 110 to 160: height 500 mm)	B29745 B29747 B2974			749		
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746 B29748 B29			750		
Subframe, adjustable in height, height 130 mm, for example for units with fresh air filter	B33657 B33659 B3			B33	661	
FDA conforming software AtmoCONTROL (FDA edition). Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit. Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1					
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence	FDAQ2					
DAkkS calibration for one free-selectable temperature value according to method C (DKD-R 5-7)	E39696					
DAkkS calibration for further temperature values according to method C (DKD-R 5-7)	E39697					
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer	D00124					
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 9 measuring points (size 30), 27 measuring points (sizes 55 - 1060) to DIN 12880:2007-05. PQ check list as support for validation by customer. 305 € for further temperature values	D00125 D00127					
On-site IQ/OQ for a freely selectable temperature value, including temperature distribution survey for 9 measuring points (size 30), 27 measuring points (sizes 55 - 160) to DIN 12880: 2007-05 (excluding travel costs, not subject to discount, GER, AT, FR only)	DLQ100					

Accessories	30	55	75	110	160			
Extension of DLQ100 by an additional freely selectable temperature value (not subject to discount)	DLQ100A							
Individual on-site Performance Qualification (PQ)			DLQ200					
Maintenance UIS - carrying out and documentation according to Memmert maintenance plan (excluding travel costs, not subject to discount, GER, AT, FR only)	S00311							
Maintenance contract UIS - carrying out and documentation according to Memmert maintenance plan, minimum duration 3 years (excluding travel costs, not subject to discount, GER, AT, FR only)	3 S00311J							
Calibration of one freely selectable temperature value (excluding travel costs, not subject to discount, GER, AT, FR only) S00205								
Calibration of an additional temperature value (not subject to discount)			S00215					